

The StimRouter system delivers stimulation from an external pulse transmitter, or EPT, through the skin to the implanted StimRouter lead. The EPT is powered by a rechargeable lithium polymer battery. Proper care of the EPT battery will extend the battery life and ensure that the StimRouter system is able to deliver effective stimulation. Please follow the steps below, as well as any other guidance and best practices suggestions from the StimRouter team.

- **Please charge the EPT battery when it is not in use. Charging will ensure that the battery is never completely drained.**
- **Please take 15 minute breaks from stimulation every 4 hours and completely remove the electrode and EPT from the skin. Be sure to charge the EPT during breaks from stimulation.**
- **Please familiarize yourself with battery indicators.**



Yellow flashing Patient Programmer icon on the Patient Programmer:

Patient Programmer battery is low – charge as soon as possible.



While Charging, Patient Programmer Icon GLOWS RED and “E” Appears in the Digital Display:

Charging Error - Reconnect the charger cable to the Patient Programmer, or follow instructions in the Maintenance and Cleaning chapter of the User’s Guide for replacing the Patient Programmer Battery.



EPT icon flashes yellow:

EPT battery is low – charge as soon as possible.

- **Before placing the electrode or beginning stimulation, please ensure proper skin preparation steps are followed. Improper skin care can cause the StimRouter System to work harder and drain the EPT battery faster.**
 - Please refer to the StimRouter User Guide for instructions for proper skin preparation.
 - Place the electrode on the skin.
 - Wait 2 minutes after placing the electrode on the skin before starting stimulation. The gel pads of the electrode use water to conduct electrical stimulation. Waiting for the gel pads to absorb some of the water on the skin will rehydrate the pads and ensure a better connection is made between the skin and electrode.

More information about StimRouter can be found by:

- Contacting your StimRouter representative
- Consulting the StimRouter User’s Guide
- Calling 800-211-9136
- Visiting www.StimRouter.com
- Visiting the StimRouter YouTube page

The StimRouter Neuromodulation System is indicated for pain management in adults who have severe intractable chronic pain of peripheral nerve origin, as an adjunct to other modes of therapy (e.g., medications). The StimRouter is not intended to treat pain in the craniofacial region.

Do not use the StimRouter Neuromodulation System in users who have an implanted demand cardiac pacemaker, implanted cardioverter defibrillator (ICD), or other implanted active device, or who have bleeding disorders that cannot be stopped in advance of the StimRouter implantation procedure. Do not use the system where a metallic implant or a cancerous lesion is present in the immediate implant area. Effects of stimulation during pregnancy are not known. StimRouter is capable of producing skin irritation and muscle ache in the area of stimulation.

Full prescribing information can be found in the Clinician Guide or <https://stimrouter.com/safety-information>.

Bioventus and the Bioventus logo are registered trademarks of Bioventus LLC.

StimRouter is a registered trademark of Bioness Inc. | www.stimrouter.com | **Rx Only**

©2022 Bioventus LLC

